



## INFORMATION DATA SHEET

### Anti-Endomysial Antibody Test System

Catalogue No: K4818E, K9618E, K5018E, K0018E, K25018E

Containing the following kit components:

1. Monkey Oesophagus (Endomysial) slides	S4205E, S8205E, S5205E or S0205E
2. Universal negative control	C000N/C000N-0.5
3. Endomysial IgA positive control	C026N/CO26N-0.5/C026S
4. FITC IgA Anti-human Conjugate	J527
5. PBS Buffer	R002
6. Mounting medium	R005

Each kit component has been assessed and has been found to be non-hazardous.

However it should be noted that controls contain human sera and although they have been screened and found to be negative for HBsAg and antibodies to HIV-1, HIV-2, P24 antigen and HCV, these tests cannot guarantee the absence of infectious agents. All human components should be handled with appropriate care. The controls also contain 0.01% Thiomersal as a preservative. Although this is at a low concentration, these reagents should be considered toxic. They should not be ingested or allowed to come into contact with either the skin or the mucous membranes.

The information in this data sheet is to our best knowledge correct and complete and is offered in good faith as accurate. It is the user's obligation to evaluate and use this product safely and to comply with all applicable laws and regulations.

Bio-Diagnostics Ltd  
Upton Industrial Estate  
Rectory Road  
Upton upon Severn  
Worcestershire  
WR8 0XL  
United Kingdom

Telephone : +44 (0)1684 592262

Email: enquiries@bio-diagnostics.co.uk

Date:

14/03/2017