



INFORMATION DATA SHEET

Anti-Ovary/Testes/Adrenal Antibody Test System

Catalogue No: K4006

Containing the following kit components:

1. Monkey Ovary/Testes/Adrenal slides	T4006/T8006
2. Universal negative control	C000N/C000N-0.5
3. Ovary (theca cells) antibody Positive Control no:	C030N/C030N-0.5
4. Adrenal antibody Positive Control no:	C017N/C017N-0.5
5. FITC primate with EB	J502
6. PBS Buffer	R002
7. Mounting medium	R005

Each kit component has been assessed and has been found to be non-hazardous.

However, it should be noted that controls contain human sera and although these have been tested for Hepatitis B and HIV, these components should be treated as potentially infectious. The conjugate and control included in the kit contain 0.1% sodium azide, 0.01% thiomersal or 0.05% ProClin as preservatives. Although these are at low concentrations, these reagents should be considered toxic. They should not be ingested or allowed to come into contact with either the skin or the mucous membranes.

The information in this data sheet is to our best knowledge correct and complete and is offered in good faith as accurate. It is the user's obligation to evaluate and use this product safely and to comply with all applicable laws and regulations. Please also refer to the Warnings and Precautions section in the Instructions for Use supplied in the kit.

Bio-Diagnostics Ltd
Upton Industrial Estate
Rectory Road
Upton upon Severn
Worcestershire
WR8 0XL
United Kingdom
Telephone : +44 (0)1684 592262
Email: enquiries@bio-diagnostics.co.uk

Date: 05/02/2020