



## INFORMATION DATA SHEET

### Anti-Adrenal Antibody Test System

Catalogue No: K4813 & K9613

Containing the following kit components:

- |                               |                         |
|-------------------------------|-------------------------|
| 1. Monkey Adrenal slides      | S4209/S8209/S5209/S0209 |
| 2. Universal negative control | C000N/C000N-0.5         |
| 3. Adrenal positive control   | C017N/C017N-0.5         |
| 4. FITC primate with EB       | J502                    |
| 5. PBS Buffer                 | R002                    |
| 6. Mounting medium            | R005                    |
| 7. Blotters                   |                         |
| 8. Coverslips                 |                         |

Each kit component has been assessed and has been found to be non-hazardous.

All human components have been tested for (HBsAg) and HTLVIII/LAV and found to be negative. However, this does not assure the absence of (HBsAg) or HTLVIII/LAV. All human components should be handled with appropriate care.

The controls included in the kit contain 0.01% thiomersal as preservatives. Although these are at low concentrations, these reagents should be considered toxic. They should not be ingested or allowed to come into contact with either the skin or the mucous membranes.

The information in this data sheet is to our best knowledge correct and complete and is offered in good faith as accurate. It is the user's obligation to evaluate and use this product safely and to comply with all applicable laws and regulations. Please also refer to the Warnings and Precautions section in the Instructions for Use supplied in the kit.

Bio-Diagnostics Ltd  
Upton Industrial Estate  
Rectory Road  
Upton upon Severn  
Worcestershire  
WR8 0XL  
United Kingdom  
Telephone : +44 (0)1684 592262  
Email: enquiries@bio-diagnostics.co.uk

Date: 02/10/2019