



## **INFORMATION DATA SHEET**

### **ANA Centromere Positive Control**

Catalogue No: C018N; C018N-0.5; C018S

These preparations have been assessed and have been found to be non-hazardous.

This positive control has been tested and found negative for STS, HIV1 Ag, HIV1/2 antibody, HCV antibody and HCV PCR (NAT) by FDA approved tests. However, these tests cannot guarantee the absence of infectious agents. All human components should be handled with appropriate care.

The control contains 0.01% thiomersal as a preservative. Although this is at a low concentration, these reagents should be considered toxic.

The information in this data sheet is to our best knowledge correct and complete and is offered in good faith as accurate. It is the user's obligation to evaluate and use this product safely and to comply with all applicable laws and regulations.

Bio-Diagnostics Ltd  
Upton Industrial Estate  
Rectory Road  
Upton upon Severn  
Worcestershire  
WR8 0XL  
United Kingdom  
Telephone : +44 (0)1684 592262  
Email: enquiries@bio-diagnostics.co.uk

Date: 12 September 2023